

Mini gel filtration chromatography for the MCMs

Preparation of Sephacryl 300 HR

1. Bring S300 resin to room temp.
2. Gently resuspend and pour ~2 mL into a 10 mL plastic column. Drain off ethanol.
3. Add 2 mL room temp S/0.15 buffer, resuspend gently, drain. Repeat several times.
4. Resuspend resin to ~50% in S/0.15

Preparation of 1 mL pipette column

1. Trim a tiny bit of the tip of a 1 mL glass pipet so that a syringe cap fits. Trim the top so a round gel loading tip can reach the 0 mL mark.
2. Insert a small plug of silanized glass wool into the top of the pipette and force into tip with compressed air.
3. Place a piece of tubing (2 inches) on top to act as a reservoir.

Pouring the column

1. Run 1 mL S/0.15 through the pipette column, cap the bottom, and fill with S/0.15 to the reservoir.
2. Add 50% S300 to the reservoir and let slowly settle into the pipette. Add more resin as needed to bring the column volume up to ~1 mL.
3. Let several bed volumes of S/0.15 drip through the column to help pack the resin.
4. Always keep the column topped off with buffer when not in use, cap the bottom and use a binder clip to seal the reservoir on top.

Calibration

1. To find the void volume, dissolve a small amount of blue dextran in S/0.15, load 25 μ L and let it enter the column. Chase with buffer to determine how much is needed to elute the blue dex peak (i.e. the void volume).
2. Next, load a 40 μ L cocktail of protein standards (10 μ L each of 8 mg/mL thyroglobulin, 10 mg/mL apoferritin, 1 mg/mL β -amylase, and 1 mg/mL BSA) and let it enter the column. Chase with buffer and collect one-drop (~21 μ L/drop) fractions starting near the void volume.
3. Add 4 μ L 6x loading dye to fractions, boil, and run 10 μ L of each on one or two 12% SDS-PAGE gels at 200V.
4. Sypro stain and quantitate the bands to find the peak fractions of each.

Running the MCM samples

1. Run 1-2 bed volumes of S/0.15 through the column to clean it out.
2. Load 20-25 μ L purified MCM sample and let it enter the column.
3. Add a bed volume of buffer and collect 20 one-drop fractions beginning near the void.
4. Add 4 μ L 6x loading dye to each fraction, boil, and run 10 μ L of each on 7% SDS-PAGE gels at 200V.
5. Sypro stain and quantitate the bands to find the peak(s).

MCM2-7 and 467 hexamers run slightly larger than thyroglobulin. MCM3 runs as a dimer.

Gel filtration chromatography

A 1 mL glass column was packed with Sephacryl 300 HR (Sigma) equilibrated in buffer S/0.15 and calibrated with a standard molecular weight markers including blue dextran (2,000 kDa), thyroglobulin (669 kDa), apoferritin (443 kDa), β -amylase (200 kDa), and BSA (66 kDa). Small samples of purified protein (20-25 μ L) were subjected to analytical gel filtration chromatography run by gravity flow at room temperature. 21 μ L fractions were collected and analyzed by SDS-PAGE and staining with Sypro orange.